

sightsense

Entry point

1Home — household overviewboth modes/home

What this page shows

This is the landing screen after login. Shows the whole household at a glance — current mode (night/day), active alerts, and system health.

Resident A

Mode indicatorLast sleep scoreCurrent room zoneActive alerts

Individual overview modes/resident/:id

Components

Current mode (night training / day detection)Attribution confidence todayLast night sleep scoreBiometric profile status (enrolled / calibrating)Weekly trend summaryActive alerts for this personCurrent room zoneActivity state (static / dynamic / absent)

Entry point into all sub-pages for this resident. Navigation sidebar lets family switch between Sleep, Activity, Safety, Report, Attribution pages — all scoped to this person.

3Sleep & breathingnight mode only/resident/:id/sleep

Components

Nightly breathing rate chartSleep quality score per nightApnea event log (rate <6 brpm or gap >10s)Restlessness index8-week trend with drift markersVs personal baseline comparison

Attribution confidence per night

Strongest page technically. Overnight is your highest-confidence window. Everything here is trained data, not inferred.

4Activity & routineboth modes/resident/:id/activity

Night mode (attributed)

Routine heatmap (solo window hours) Routine regularity score vs baseline
Prolonged inactivity in own room

Day mode (detection only — no identity)

Zone activity log (which room, not who) Shared space occupancy periods
Multi-occupant windows flagged

Day mode activity shown per zone, not per person. Label it clearly: "Room A activity" not "Resident A activity."

5Live sensing viewboth modes/resident/:id/live

Components

Room floor plan with sensor positions Live CSI waveform (raw → filtered toggle)
Breathing rate live chart (night only) Four-pillar Bayesian confidence meter
Activity state badge (day) Occupancy confidence (day) Mode indicator (night training / day detection)

This is your demo centrepiece. During night mode it shows live breathing. During day mode it shows activity state and zone occupancy only — vitals section visibly greyed out with label "Day mode — training paused." That greying out is a feature to show judges, not a weakness to hide.

6Biometric profilenight mode — trained data/resident/:id/profile

Components

Breathing rate distribution (enrolled) Breathing regularity distribution
Breathing amplitude distribution

Calibration progress (days since enrolment) Profile confidence quality score
Last re-enrolment date Drift from baseline (current vs enrolled)

This page makes Pillar 3 visible and tangible. Shows judges and families what "trained on this person" actually means — not a black box, a visible distribution. Differentiates from every competitor who shows no such transparency.

Safety & alerts

7SOS & safety eventsSOS/resident/:id/safety

Real-time triggers (identity-independent)

Fall detection (room-scoped, fires in any mode)Prolonged inactivity alert (>N mins no motion)Apnea alert (night mode, breathing stops >10s)Bed exit not returned (>10 mins, night)

SOS dispatch

Family push notification (immediate)Caregiver / nurse SMS or call triggerAlert log with timestamp + typeAcknowledge / false positive flagResponse time tracker

Fall detection fires regardless of mode and regardless of which person it is — room-scoped sensor is sufficient and identity-independent. Do not claim you know which resident fell in a shared space. In a shared space you know a fall happened in that zone, not who fell.

Reporting & transparency

8Weekly reporttrained data/resident/:id/report

Components

LLM-generated summary (Claude API)Sleep quality this week vs 4-week averageRoutine changes flaggedAttribution quality this week (%)Alerts triggered and resolvedExport to PDF

Share with family / clinician

9Attribution transparencyboth modes/resident/:id/attribution

Components

Four-pillar breakdown per dayHigh / partial / excluded window ratioAttribution rate 30-day trendMulti-occupant period heatmapPlain language explanation of confidenceWhat the system is and isn't confident about

Your strongest differentiator with judges. No competitor surfaces this. Make the four pillars visible, legible, and honest.

System

10Settingsboth/settings

Both resident profilesSensor pair status per roomSOS contacts (family + nurse/caregiver)Alert thresholds per residentNight mode schedule (auto or manual)Re-enrolment triggerWhat is / isn't recorded (transparency)

11Demo / judge modedemo only/demo

Raw CSI → filtered breathing waveform toggle
Simulated apnea triggerAlert firing
sequenceNight → day mode switch live
Edinburgh dataset replaySOS dispatch
simulationFour-pillar confidence animation

Isolated at /demo. Clearly labelled. Never linked from main nav. Judges only.